

THE POWER OF INNOVATION.

Sterilizing in the laboratory. Easier, faster and safer with Systec.

Program overview





Systec Laboratory autoclaves

PERFORMANCE AND RELIABILITY

For your sterilization tasks. Easier, precise, more reproducible and validatable.

Made-to-measure for any laboratory application. Modular designed autoclaves.



Our aim is always to make steam sterilization in the laboratory safer, easier and more precise, reproducible and validatable. With over 25 years experience and constant cooperation with experts in practice, we know how to tackle and provide optimal solutions for even the most complex sterilization tasks.

We have the knowledge and experience to produce the results!

Our expertise and know-how are available for you worldwide through specialized partners.

Product expertise:

- Development
- Design
- Production of serial products
- Production of specialty products
- Application and technical advice


Additional services:

- Installation and start-up
- Special technical developments
- Tests and process development
- Individual service on-call
- Contract service
- Qualification and validation
- GMP-compliant documentation

We manufacture our horizontal bench-top autoclaves, horizontal floor-standing autoclaves and vertical floor-standing autoclaves in different performance categories as sterilization tasks in laboratories are highly variable.











For elementary laboratory applications with limited availability of options for process optimization.





For basic laboratory applications without availability of options for process optimization.

Sterilization made easier, precise, more reproducible and validatable with Systec laboratory autoclaves.







www.systec-lab.com www.systec-lab.com

Systec V–Series

FOR BETTER STERILIZATION.

Systec

Vertical floor-standing autoclaves Systec V-Series

- Highest requirements with regard to safe, accurate, reproducible and validatable sterilization processes
- Sterilization temperatures up to 140 °C, steam pressure up to 4 bar absolute pressure (optional 150 °C/5 bar absolute pressure for Systec VX/VE, 65 to 150 liters)
- Fully automatic, micro-processor controlled
- Up to 25 sterilization programs for different sterilization applications

Systec V-Series 8 sizes, 40 to 150 liters chamber capacity







3 model ranges with different performance characteristics







Systec V-55



Systec V-65



Systec V-75



Systec V-95 Systec V-100





Systec V-120 Systec V-150

- 4 -

8 models 40 to 150 l

chamber capacity

- 5 -

www.systec-lab.com www.systec-lab

Systec D–Series

FOR BETTER STERILIZATION.

6 models 45 to 200 l chamber capacity

Systec D10

Horizontal bench-top autoclaves Systec D-Series

- Highest requirements with regard to safe, accurate, reproducible and validatable sterilization processes
- Sterilization temperatures up to 140 °C, steam pressure up to 4 bar absolute pressure (optional 150 °C/5 bar absolute pressure for Systec DX/DE, 65 to 150 liters)
- Fully automatic, micro-processor controlled
- Up to 25 sterilization programs for different sterilization applications, self-programmable

Systec D-Series

6 sizes, 45 to 200 liters chamber capacity







3 model ranges with different performance characteristics



Systec D-45



Systec D-65



Systec D-90



Systec D-100



Systec D-150



Systec D-200

Systec

-7-

Systec H–Series

FOR BETTER STERILIZATION.



Shorter and more efficient process times
Thanks to the completely newly developed control
of the autoclaves.



Horizontal floor-standing autoclaves Systec H-Series

- Highest requirements with regard to safe, accurate, reproducible and validatable sterilization processes
- Sterilization temperatures up to 140 °C, steam pressure up to 4 bar absolute pressure (optional up to 150 °C/5 bar absolute pressure)
- Fully automatic, micro-processor controlled
- Internal memory capacity for storing up to 10 years of data backups
- Up to 100 sterilization programs for different sterilization applications, self-programmable

Systec H-Series 16 sizes, 65 to 1580 liters chamber capacity







65 to 1580 l chamber capacity

16 models

- 9 -

www.systec-lab.com www.syst

Systec H–Series 2D

FOR BETTER STERILIZATION.



Shorter and more efficient process times
Thanks to the completely newly developed control
of the autoclaves.



Horizontal pass-through autoclaves Systec H-Series 2D

- Highest requirements with regard to safe, accurate, reproducible and validatable sterilization processes
- For gas-tight installation in partition walls of biological safety laboratories or clean rooms
- Sterilization temperatures up to 140 °C, steam pressure up to 4 bar absolute pressure (optional up to 150 °C/5 bar absolute pressure)
- Fully automatic, micro-processor controlled
- Internal memory capacity for storing up to 10 years of data backups
- Up to 100 sterilization programs for different sterilization applications, self-programmable

Systec H-Series 2D 13 sizes, 90 to 1580 liters chamber capacity



OPTIONAL: Stainless steel casing





10 – – 11 –

www.systec-lab.com www.syste


STERILIZATION OF LIQUIDS.

Heating up

- The actual sterilization time of e.g. 15 to 20 minutes at 121 °C is only a fraction of the total time involved for an autoclave procedure. Especially in the case of sterilizing liquids, the heating up and cooling down phases are considerably longer.
- The conventional procedure

In previously used conventional systems, even if the intended sterilization temperature has been reached within the autoclave, the liquids to be sterilized are often only at about 100 °C; the temperature equilibrium time between chamber and liquids normally takes much longer.

■ Up to 50 % shorter heat-up times as standard Due to the combined temperature and pressure regulation, the chamber pressure is increased during the heat-up phase. The result: more rapid temperature equilibrium in the liquids and a shorter heat-up time.



Systec D-Series/H-Series

Systec V-Series

Cooling

- The cooling process for liquids is also very slow; this is because, without active rapid cooling, the heat can only be reduced to below 100 °C by dissipating the heat via the chamber insulation by radiation.
- New system and process technology now make it possible to substantially reduce the overall time required for the sterilization process. This means that several hours of time can be saved! It also means that the media is not exposed to heat unnecessarily long time.
- Systec offers many functions for its autoclaves guaranteeing safe liquid sterilization processes at higher productivity. Many of these functions are standard or available as options depending on the model range selected.

Standard functions in all models

- Temperature- and pressure-dependent door locking in line with international standards and regulations.
- Redundant process control; temperature and pressure are continuously monitored and controlled during the entire sterilization cycle.
- Rapid heat-up via optimized heat transfer to the liquid media.
- Flexible PT-100 temperature sensor for temperature measurement in a reference vessel:
 - Guarantees attainment of the desired sterilization temperature in the liquid media.

www.systec-lab.com

- Guarantees cooling of the liquid media to a temperature that is safe for removal.

- 12 - - 13 -

Systec Additional options

QUICK-COOLING OPTIONS FOR STERILIZATION OF LIQUIDS.





STERILIZATION OF SOLIDS AND WASTE IN DISPOSAL BAGS.















■ Water Cooling with support pressure provide a coil wrapped around the outside of the sterilization vessel. During the cooling phase, cool water circulates through the coils, reducing cool down time. Support pressure prevents liquids from boiling over during the cooling process, eliminating evaporative loss of sterilized liquids. Reduction of cooling time up to 60 %.






Systec D-Series/H-Series

Systec V-Series

■ Radial ventilator for air circulation and accelerated heat removal from the sterilizer chamber. Combined with cooling coils and support pressure, the radial ventilator creates a more rapid exchange of heat from the chamber by utilizing an internal fan to circulate air. Reduction of cooling time up to 70 %.





Systec D-Series/H-Series

Systec V-Series

■ Ultracooler is an additional heat exchanger behind the radial ventilator. This special ridged coil provides a substantial increase in surface area for fast and efficient removal of heat during the cooling process. Reduction of cooling time up to 90 %.

 Vacuum system removes the air highly effectively from solids, tubing, porous materials, fabrics and disposal bags; in this way, the steam is able to penetrate completely. The process includes a fractionated pre-vacuum phase in combination with the standard steam generator.







Systec D-Series/H-Series

Systec V-Series

• Superdry increases the drying efficiency for solids and porous materials such as filters and fabrics. Heat energy from the standard steam generator is transferred to the heating coils around the body of the sterilization chamber and is used for drying. Deep-vacuum drying using the optional vacuum device in conjuction with Superdry avoids the necessity for subsequent drying in a separate drying cabinet.





Systec D-Series/H-Series

Systec V-Series

■ Exhaust air filtration system for the sterilization of hazard biological substances. The filter is automatically sterilized during each sterilization process, monitored by a PT-100 temperature sensor. Through air exhaust filtration and condensate inactivation, it is ensured that no microorganisms can escape. For waste bags

- 15 -

- 14 -

Systec Mediaprep

FOR BETTER MEDIA PREPARATION AND STERILIZATION.

Media preparators Systec Mediaprep

- Gentle preparation, rapid cooling
- Continuous mixing
- Sterilization temperature up to 136 °C
- Steam pressure up to 3,5 bar absolute pressure



Systec Mediaprep 7 sizes, 10 to 120 liters chamber capacity



7 models 10 to 120 l chamber capacity

www.systec-lab.com www.systec-lab.

Systec Media fill

FOR AUTOMATED POURING AND STACKING OF CULTURE MEDIA.



Plate pourer and tube filler Systec Mediafill





- For precise, automated and unattended pouring of culture media into Petri dishes, bi-plates, tri-plates and test tubes
- Touch display for easy and intuitive programming
- Even surfaces, no holes or gaps, easy to clean
- Three stacker carousel heights for up to 220, 440 and 660 Petri dishes



8 -



Systec GmbH
Konrad-Adenauer-Straße 15
35440 Linden, Germany
T +49 6403 670 70 0
F +49 6403 670 70 222
info@systec-lab.de
www.systec-lab.com



Subsidiary Switzerland:
Systec Schweiz GmbH
Gewerbestrasse 8
CH-6330 Cham, Switzerland
T +41 41 781 52 80
F +41 41 781 52 79
info@systec-lab.ch
www.systec-lab.ch



Subsidiary China:
Systec (Shanghai) trading Co., Ltd
D1-201, No.6000 Shenzhuan Rd,
Songjiang 201619
Shanghai, China
T +86 21 6019 0256
info@systec-lab.com.cn
www.systec-lab.com.cn

Scan to visit our website:

